



LG

Life's Good


MEDICAL MONITOR SERVICE MANUAL

CHASSIS : LM2HA

MODEL : 55MH5K 55MH5K-WX

CAUTION

BEFORE SERVICING THE CHASSIS, READ THE SAFETY PRECAUTIONS IN THIS MANUAL.



P/NO : MFL71804643 (2304-REV00)

Any reproduction, duplication, distribution (including by way of email, facsimile or other electronic means), publication, modification, copying or transmission of this Service Manual is STRICTLY PROHIBITED unless you have obtained the prior written consent of the LG Electronics entity from which you received this Service Manual. The material covered by this prohibition includes, without limitation, any text, graphics or logos in this Service Manual.

CONTENTS

| | |
|-------------------------------------|------------------|
| CONTENTS | 2 |
| SAFETY PRECAUTIONS | 3 |
| SPECIFICATION | 4 |
| BLOCK DIAGRAM | 10 |
| EXPLODED VIEW | 11 |
| TROUBLE SHOOTING GUIDE | 12 |
| DISASSEMBLE GUIDE | APPENDIX1 |
| FIRMWARE UPDATE GUIDE | APPENDIX2 |

SAFETY PRECAUTIONS

IMPORTANT SAFETY NOTICE

Many electrical and mechanical parts in this chassis have special safety-related characteristics. These parts are identified by \triangle in the Exploded View.

It is essential that these special safety parts should be replaced with the same components as recommended in this manual to prevent Shock, Fire, or other Hazards.

Do not modify the original design without permission of manufacturer.

General Guidance

An **isolation Transformer should always be used** during the servicing of a receiver whose chassis is not isolated from the AC power line. Use a transformer of adequate power rating as this protects the technician from accidents resulting in personal injury from electrical shocks.

It will also protect the receiver and its components from being damaged by accidental shorts of the circuitry that may be inadvertently introduced during the service operation.

If any fuse (or Fusible Resistor) in this TV receiver is blown, replace it with the specified.

When replacing a high wattage resistor (Oxide Metal Film Resistor, over 1 W), keep the resistor 10 mm away from PCB.

Keep wires away from high voltage or high temperature parts.

Before returning the receiver to the customer,

always perform an **AC leakage current check** on the exposed metallic parts of the cabinet, such as antennas, terminals, etc., to be sure the set is safe to operate without damage of electrical shock.

Leakage Current Cold Check(Antenna Cold Check)

With the instrument AC plug removed from AC source, connect an electrical jumper across the two AC plug prongs. Place the AC switch in the on position, connect one lead of ohm-meter to the AC plug prongs tied together and touch other ohm-meter lead in turn to each exposed metallic parts such as antenna terminals, phone jacks, etc.

If the exposed metallic part has a return path to the chassis, the measured resistance should be between 1 M Ω and 5.2 M Ω .

When the exposed metal has no return path to the chassis the reading must be infinite.

An other abnormality exists that must be corrected before the receiver is returned to the customer.

Leakage Current Hot Check (See below Figure)

Plug the AC cord directly into the AC outlet.

Do not use a line Isolation Transformer during this check.


Connect 1.5 K / 10 watt resistor in parallel with a 0.15 uF capacitor between a known good earth ground (Water Pipe, Conduit, etc.) and the exposed metallic parts.

Measure the AC voltage across the resistor using AC voltmeter with 1000 ohms/volt or more sensitivity.

Reverse plug the AC cord into the AC outlet and repeat AC voltage measurements for each exposed metallic part. Any voltage measured must not exceed 0.75 volt RMS which corresponds to 0.5 mA.

In case any measurement is out of the limits specified, there is possibility of shock hazard and the set must be checked and repaired before it is returned to the customer.

Leakage Current Hot Check circuit



When 25A is impressed between Earth and 2nd Ground for 1 second, Resistance must be less than 0.1

*Base on Adjustment standard

SPECIFICATION

NOTE : Specifications and others are subject to change without notice for improvement.

1. General Specification

| No. | Item | | Content | | Remark |
|------|-----------------|--------------------|--|--|--|
| 1 | Customer | | BRAND | | BRAND : LGE Product |
| 2 | User Model Name | | 55MH5K | | |
| 3 | Sale region | | Worldwide | | |
| 4 | Feature | | 55" UHD LCD Medical Monitor (with Direct Bonding Glass) | | |
| 5 | Chassis Name | | LM2HA | | "PCB, Product spec., Adjustment spec." |
| 6 | General Scope | External SW & Adj. | Input, PIP/PBP, Hotkey1, Hotkey2, Brightness-, Brightness+, OK, Contrast-, Contrast+, Menu | | 10 Tack Key |
| | | Function | Picture Mode, Ratio, S.E.S, Six Color, LG Calibration Studio, Factory calibration, PBP/PIP, RS-232C, User Preset | | |
| 7 | Input / Output | SDI IN | 4ea (12G-SDI) | | |
| | | SDI OUT | 4ea (12G-SDI) | | |
| | | SFP | Support 1Rx SFP module | | Not supported 2Rx & TxRx SFP module |
| | | DVI-D IN | 1ea | | |
| | | DVI-D OUT | 1ea | | |
| | | HDMI IN | 1ea (HDMI2.0, UHD, Support HDR) | | |
| | | DP IN | 1ea (DP1.2, UHD, Not support HDR) | | |
| | | DP OUT | 1ea (DP1.2, UHD, Not support HDR) | | |
| | | USB UP | 1ea (USB3.0 Upstream) | | |
| | | USB IN | 1ea (USB3.0 Downstream) | | USB downstream Port(Up to 5V 0.9A DC) |
| | | SERVICE ONLY | 1ea | | |
| | | RS232C | 1ea | | |
| | | AC IN | 1ea | | |
| POAG | 1ea | | | | |
| | | SOUND | N/A (No speaker, No HP out) | | |
| 8 | Power Code | | Length | 2.8m ± 0.05 | Refer to Suffix standard and power cord Table. |
| | | | Weight | 378g ± 10% | |
| | | | Color | White | |
| 9 | Cable | DVI | Length | 1.8m ± 0.05m | Not support |
| | | | Weight | 0.11kg ± 10%, | |
| | | | Color | Black | |
| | | | Impedance | CABLE(Differential) Impedance : 100Ω±10Ω | |
| | | DisplayPort | Length | 1.80m±0.05m | |
| | | | Weight | 0.10kg ± 10% | |
| | | | Color | Black | |
| | | | Impedance | CABLE(Differential) Impedance : 100Ω±10Ω | |

| No. | Item | | Content | | Remark |
|-----|---------------|------|---|---------------|--------|
| 9 | Cable | HDMI | Length | 1.80m ± 0.05m | |
| | | | Weight | 0.09kg ± 10% | |
| | | | Color | Black | |
| | | | CABLE(Differential) Impedance: 100Ω±10Ω | | |
| 10 | Power | PSU | Input: 100-240 V~ 50-60 Hz 2.5 – 1.25A Output: 20.5V, 9.0A | | |
| 11 | Power Rating | | Input : AC100~240V 50~60Hz, 24V 7.5A | | |
| 12 | Accessory Box | | Weight : 0.8kg | | |

2. Engineering Specification

| No. | Item | Spec | | Remarks | |
|-----|---------------------|---|------------|------------------------------|--|
| 1 | Operating Frequency | SDI | Horizontal | N/A | None EDID. Refer to Support resolution |
| | | | Vertical | N/A | |
| | | DVI | Horizontal | 30 - 83 kHz | |
| | | | Vertical | 56 - 61 Hz | |
| | | HDMI | Horizontal | 30 - 135 kHz | |
| | | | Vertical | 56 - 61 Hz | |
| | | DP | Horizontal | 30 - 135 kHz | |
| | | | Vertical | 56 - 61 Hz | |
| 2 | Resolution | SDI | Max. | 12G-SDI (4096x2160 @60Hz) | 12G-SDI : 3840x2160 @60Hz (Recommend) 3G-SDI : 1920x1080 @60Hz |
| | | | Recommend | 12G-SDI (3840x2160 @60Hz) | |
| | | DVI | Max. | 1920x1080 @ 60Hz | |
| | | | Recommend | 1920x1080 @ 60Hz | |
| | | HDMI | Max. | 3840×2160 @ 60Hz | |
| | | | Recommend | 3840×2160 @ 60Hz | |
| | | DP | Max. | 3840×2160 @ 60Hz | |
| | | | Recommend | 3840×2160 @ 60Hz | |
| 3 | Input Voltage | Voltage: AC100~240V 50~60Hz, 20.5V 9.0A | | PSU | |

3. Signal Timing (Resolution)

3.1. Signal (Video & Sync)

3.1.1. H/V Timing : DisplayPort

| No. | Section | Pol. | Dot Clock [MHz] | Frequency [kHz]/[Hz] | Total Cycle (E) | Display (A) | Front Porch(B) | Sync. (D) | Back Porch(F) | Resolution | Remark |
|------------------|-----------|------|-----------------|----------------------|-----------------|-------------|----------------|-----------|---------------|-------------|--------|
| 1 | H(Pixels) | - | 25.175 | 31.469 | 800 | 640 | 16 | 96 | 48 | 640 x 480 | |
| | V(Lines) | - | | 59.94 | 525 | 480 | 10 | 2 | 33 | | |
| 2 | H(Pixels) | + | 40 | 37.879 | 1056 | 800 | 40 | 128 | 88 | 800 x 600 | |
| | V(Lines) | + | | 60.317 | 628 | 600 | 1 | 4 | 23 | | |
| 3 | H(Pixels) | - | 65 | 48.363 | 1344 | 1024 | 24 | 136 | 160 | 1024 x 768 | |
| | V(Lines) | - | | 60 | 806 | 768 | 3 | 6 | 29 | | |
| 4 | H(Pixels) | + | 79.99 | 54.347 | 1472 | 1152 | 32 | 96 | 192 | 1152 x 864 | |
| | V(Lines) | + | | 60.05 | 905 | 864 | 1 | 3 | 37 | | |
| 5 | H(Pixels) | + | 74.25 | 45 | 1650 | 1280 | 110 | 40 | 220 | 1280 x 720 | |
| | V(Lines) | + | | 60 | 750 | 720 | 5 | 5 | 20 | | |
| 6 | H(Pixels) | + | 108 | 63.981 | 1688 | 1280 | 48 | 112 | 248 | 1280 x 1024 | |
| | V(Lines) | + | | 60.02 | 1066 | 1024 | 1 | 3 | 38 | | |
| 7 | H(Pixels) | + | 108 | 60 | 1800 | 1600 | 24 | 80 | 96 | 1600 x 900 | |
| | V(Lines) | + | | 60 | 1000 | 900 | 1 | 3 | 96 | | |
| 8 | H(Pixels) | + | 148.5 | 67.5 | 2200 | 1920 | 88 | 44 | 148 | 1920 x 1080 | |
| | V(Lines) | - | | 60 | 1125 | 1080 | 4 | 5 | 36 | | |
| 9 | H(Pixels) | + | 241.5 | 88.79 | 2720 | 2560 | 48 | 32 | 80 | 2560 x 1440 | |
| | V(Lines) | - | | 59.95 | 1481 | 1440 | 3 | 5 | 33 | | |
| 10 | H(Pixels) | + | 266.64 | 66.66 | 4000 | 3840 | 8 | 144 | 8 | 3840 x 2160 | |
| | V(Lines) | - | | 30 | 2222 | 2160 | 54 | 5 | 3 | | |
| 11 ¹⁾ | H(Pixels) | + | 533.28 | 133.32 | 4000 | 3840 | 48 | 32 | 80 | 3840 x 2160 | |
| | V(Lines) | - | | 60 | 2222 | 2160 | 54 | 5 | 3 | | |

* 1) : DP Preferred timing

3.1.2. H/V Timing : HDMI / DVI

| No. | Section | Pol. | Dot Clock [MHz] | Frequen- cy [kHz]/[Hz] | Total Cycle (E) | Dis- play (A) | Front Porch(B) | Sync. (D) | Back Porch(F) | Resolution | Remark |
|------------------|-----------|------|-----------------|---------------------------|-----------------|------------------|----------------|-----------|---------------|-------------|--------|
| 1 | H(Pixels) | - | 25.175 | 31.469 | 800 | 640 | 16 | 96 | 48 | 640 x 480 | |
| | V(Lines) | - | | 59.94 | 525 | 480 | 10 | 2 | 33 | | |
| 2 | H(Pixels) | + | 40 | 37.879 | 1056 | 800 | 40 | 128 | 88 | 800 x 600 | |
| | V(Lines) | + | | 60.317 | 628 | 600 | 1 | 4 | 23 | | |
| 3 | H(Pixels) | - | 65 | 48.363 | 1344 | 1024 | 24 | 136 | 160 | 1024 x 768 | |
| | V(Lines) | - | | 60 | 806 | 768 | 3 | 6 | 29 | | |
| 4 | H(Pixels) | + | 79.99 | 54.347 | 1472 | 1152 | 32 | 96 | 192 | 1152 x 864 | |
| | V(Lines) | + | | 60.05 | 905 | 864 | 1 | 3 | 37 | | |
| 5 | H(Pixels) | + | 74.25 | 45 | 1650 | 1280 | 110 | 40 | 220 | 1280 x 720 | |
| | V(Lines) | + | | 60 | 750 | 720 | 5 | 5 | 20 | | |
| 6 | H(Pixels) | + | 108 | 63.981 | 1688 | 1280 | 48 | 112 | 248 | 1280 x 1024 | |
| | V(Lines) | + | | 60.02 | 1066 | 1024 | 1 | 3 | 38 | | |
| 7 | H(Pixels) | + | 108 | 60 | 1800 | 1600 | 24 | 80 | 96 | 1600 x 900 | |
| | V(Lines) | + | | 60 | 1000 | 900 | 1 | 3 | 96 | | |
| 8 ¹⁾ | H(Pixels) | + | 148.5 | 67.5 | 2200 | 1920 | 88 | 44 | 148 | 1920 x 1080 | |
| | V(Lines) | - | | 60 | 1125 | 1080 | 4 | 5 | 36 | | |
| 9 | H(Pixels) | + | 241.5 | 88.79 | 2720 | 2560 | 48 | 32 | 80 | 2560 x 1440 | |
| | V(Lines) | - | | 59.95 | 1481 | 1440 | 3 | 5 | 33 | | |
| 10 | H(Pixels) | + | 297 | 67.5 | 4400 | 3840 | 176 | 88 | 296 | 3840 x 2160 | |
| | V(Lines) | - | | 30 | 2250 | 2160 | 8 | 10 | 72 | | |
| 11 ²⁾ | H(Pixels) | + | 594 | 135 | 4400 | 3840 | 176 | 88 | 296 | 3840 x 2160 | |
| | V(Lines) | - | | 60 | 2250 | 2160 | 8 | 10 | 72 | | |

* (1) : DVI Preferred timing.

* (2) : HDMI Preferred timing.

* Mode 9,10,11 can be supported to only HDMI.

3.1.3. Video input: DP / HDMI / DVI


| EDID 2Page Video Short Block Description | | | | | | | | | | |
|--|-----------------|---------------|-----|------------|-------|------------|-----------------|---------------------|---------------------|------------|
| EDID NO. | Horizontal size | Vertical size | P/I | V-freq(Hz) | Ratio | DVI (FHD) | HDMI1.4 (UHD30) | HDMI2.0 (UHD60 HDR) | HDMI2.0 (UHD60 SDR) | DP (UHD60) |
| 1 | 640 | 480 | P | 59.94/60 | 4:03 | | ■ | ■ | ■ | ■ |
| 2 | 720 | 480 | P | 59.94/60 | 4:03 | | | ■ | ■ | |
| 3 | 720 | 480 | P | 59.94/60 | 16:09 | ■ | ■ | ■ | ■ | |
| 4 | 1280 | 720 | P | 59.94/60 | 16:09 | ■ | ■ | ■ | ■ | ■ |
| 5 | 1920 | 1080 | I | 59.94/60 | 16:09 | ■ | ■ | ■ | ■ | |
| 90(10) | 1920 | 1080 | P | 59.94/60 | 16:09 | ■ (Native) | ■ (Native) | ■ (Native) | ■ (Native) | ■ (Native) |
| 12 | 720 | 576 | P | 50 | 16:09 | ■ | ■ | ■ | ■ | |
| 13 | 1280 | 720 | P | 50 | 16:09 | ■ | ■ | ■ | ■ | |
| 14 | 1920 | 1080 | I | 50 | 16:09 | ■ | ■ | ■ | ■ | |
| 1F | 1920 | 1080 | P | 50 | 16:09 | ■ | ■ | ■ | ■ | |
| 20 | 1920 | 1080 | P | 23.97/24 | 16:09 | | | ■ | ■ | |
| 22 | 1920 | 1080 | P | 29.97/30 | 16:09 | | ■ | ■ | ■ | |
| 5D | 3840 | 2160 | P | 23.98/24 | 16:09 | | | ■ | ■ | |
| 5E | 3840 | 2160 | P | 25 | 16:09 | | | ■ | ■ | |
| 5F | 3840 | 2160 | P | 29.97/30 | 16:09 | | | ■ | ■ | |
| 60 | 3840 | 2160 | P | 50 | 16:09 | | | ■ | ■ | |
| 61 | 3840 | 2160 | P | 59.94/60 | 16:09 | | | ■ | ■ | |

3.1.4. SDI Input(with SFP) for 32HL714S (max 12G-SDI)

- Equipment: Omnitek 4K Tool Box
- Supported Timing : max 12G-SDI

| SMPTE | Interface | Signal Type | Image Format | Signal Format | Sampling Structure | Pixel Depth | Frame Rate | SDI | SFP | Remark |
|---------------------|---------------------|----------------------|--------------|------------------------------------|-------------------------------|-------------------------------|------------------------------------|-----|-----|--|
| 259M | SD-SDI | | 720 x 486 1) | 4:2:2 | YCbCr | 10bit | 59.94i | O | O | |
| | | | 720 x 576 | 4:2:2 | YCbCr | 10bit | 50i | O | O | |
| 292M | HD-SDI | | 228 x 1080 | 4:2:2 | YCbCr | 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| | | | | | | | 60i, 59.94i, 50i, 48i, 47.95i | O | O | |
| | | | 1920 x 1080 | 4:2:2 | YCbCr | 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| | | | | | | | 60i, 59.94i, 50i, 48i, 47.95i | O | O | |
| 1280 x 720 | 4:2:2 | YCbCr | 10bit | 60p, 59.94p, 50p, 30p, 29.97p, 25p | O | O | | | | |
| ST 425-1 | 3G-SDI | Level A | 1280 x 720 | 4:2:2 | RGB / YCbCr | 10bit | 60p, 59.94p, 50p, 30p, 29.97p, 25p | O | O | |
| | | | | | | | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| | | Level A / Level B-DL | 228 x 1080 | 4:2:2 | RGB / YCbCr | 12bit / 10bit | 60i, 59.94i, 50i, 48i, 47.95i | O | O | |
| | | | | | | | 60p, 59.94p, 50p, 48p, 47.95p | O | O | |
| | | | 1920 x 1080 | 4:2:2 | RGB / YCbCr | 12bit / 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| | | | | | | | 60i, 59.94i, 50i, 48i, 47.95i | O | O | |
| | | 4:2:2 | YCbCr | 10bit | 60p, 59.94p, 50p, 48p, 47.95p | O | O | | | |
| | | | | | 12bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | | |
| | | 60i, 59.94i, 50i | O | O | | | | | | |
| | | | O | O | | | | | | |
| ST 425-5 | Quad-Link 3G-SDI 2) | Level A / Level B-DL | 4096 x 2160 | 4:2:2 | RGB & YCbCr | 12bit / 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | | 2-sample interleave division square division |
| | | | | | | | 60p, 59.94p, 50p | O | | |
| | | | 3840 x 2160 | 4:2:2 | RGB & YCbCr | 12bit / 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | | |
| | | | | | | | 60p, 59.94p, 50p | O | | |
| | | | 12bit | 30p, 29.97p, 25p, 24p, 23.98p | O | | | | | |
| | | | | O | | | | | | |
| ST 2081-10 (Mode 1) | 6G-SDI | | 4096 x 2160 | 4:2:2 | YCbCr | 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| | | | 3840 x 2160 | 4:2:2 | YCbCr | 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| ST 2082-10 (Mode 1) | 12G-SDI | | 4096 x 2160 | 4:2:2 | RGB & YCbCr | 12bit / 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| | | | | | | | 60p, 59.94p, 50p, 48p, 47.95p | O | O | |
| | | | 3840 x 2160 | 4:2:2 | RGB & YCbCr | 12bit / 10bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | |
| | | | | | | | 60p, 59.94p, 50p, 48p, 47.95p | O | O | Recommended 60P |
| | | | 12bit | 30p, 29.97p, 25p, 24p, 23.98p | O | O | | | | |

BLOCK DIAGRAM



EXPLODED VIEW


IMPORTANT SAFETY NOTICE

Many electrical and mechanical parts in this chassis have special safety-related characteristics. These parts are identified by \triangle in the EXPLODED VIEW.

It is essential that these special safety parts should be replaced with the same components as recommended in this manual to prevent Shock, Fire, or other Hazards.

Do not modify the original design without permission of manufacturer.


* When module or mainboard is defected by any reason, monitor should be changed to new one.




TROUBLE SHOOTING GUIDE

First of all, Check whether there is SVC Bulletin in GSCS System for these model.


1. Breakdown repair : No Power




2. Breakdown repair : No Raster




3. Breakdown repair : No Raster - HDMI




4. Breakdown repair : No Raster – DVI In / Out



5. Breakdown repair : No Raster – DP In / Out




6. Breakdown repair : No Raster – SDI In/Out




DISASSEMBLE GUIDE

Monitor Exploded Method


[Notice] The image quality correction value is not reflected when replacing the module or main board separately. If the module or main board is defective, please exchange the product.




1. Disassemble the screws securing the rear door.
2. Disassemble the Back Door.



3. Disassemble the VESA screw and B/C screw.
4. Disassemble the B/C Assy.



5. Disassemble the VESA Metal fixing screw.
6. Disassemble the VESA Metal.




7. Disassemble rear shield metal fixing screws, separate cables, and remove tapes.
8. Disassemble the Rear shield Metal.

Tool Description

General drivers



⚠ Caution




When removing the rear shield, lift the bottom slightly and disconnect the LVDS cable first.

Monitor Exploded Method



9. Remove the hex screws and SDI JACK fixing screws.



10. Remove the PCB screws.




▣ Tool Description

General drivers



* Indicator board rework guide



1. Remove the screw 3EA fastened to the Control PCB and remove the Control PCB.



2. Using a tool, Slightly push each key on the button (in the direction of the red arrow) so that it falls off the double-sided tape on the sheet. (Be careful not to break the button)



3. Detach the control button.



4. Detach the indicator board using the tool.

FIRMWARE UPDATE GUIDE

1. How to Monitor Firmware update by United JIG

There are two ways to update the monitor firmware.


1. How to use the update jig
2. How to use a USB memory stick

Both methods are the same, so you can choose the method that is convenient for you.

1. How to Monitor Firmware update by United JIG

- Need Software to download ISP file. And Method of JIG connection.
 - : DebugTool Driver (JIG Driver Program)
 - : ISP_Tool_V4.5.x.x.x (ISP Execute TOOL Program)
 - : SW.bin file



Step1. Connect United JIG Jig(Photo #1.) with Monitor & PC.



1. How to Monitor Firmware update by United JIG


Step2. Install MSTAR Driver

- 1) If connect USB Cable to Jig, you can see pop-up screen as below.
- 2) Find drive file of your PC Windows version and selection(install).
- 3) If finish installing driver, please re-booting PC.




1. How to Monitor Firmware update by United JIG

Step3. ISP Execution TOOL Program




1) Config (program / JIG setting)

- : Check Use "USB", "Auto release".
- : Adjust I2C Speed suitably




2) Connect (Connect JIG and Flash memory)

- : After Monitor AC Power On and connect JIG and Cable.
- : If push the Connect icon, you can see Flash memory name on the pop-up screen.
- : Push OK.



Occur error message,
re-try again.
(connection check)
(restart program)
(reboot)




1. How to Monitor Firmware update by United JIG

3) Read (prepare Bin file)

: Push "Read" icon and prepare ISP bin file. ※ Read is recommended after Connect Success.


: Check : End Addr. ※ Sometimes the value 0x1FFFF is displayed as an example, in which case you need to read it again.



4) Auto (Download)

: After check A and B and push the "Run".

: You can see download time and finish message .



2. How to Monitor Firmware update by USB Memory stick

1. Please prepare an USB memory stick. (Except for the external HDD disk)


[Note] The file system must use the following format. (NTFS, FAT, FAT32)

2. Please make a new folder like as LG_MNT of USB memory stick. (Ex. E:\WLG_MNT)

3. Please copy monitor FW(MH5K_MSTAR_MST9W00V4_V3.X.X_XXXX_revXXXX_XXXXXX.bin) to the folder

[Note] The downloading feature does not work properly when using the download file from the company security.

4. Please insert the USB memory stick to the "SERVICE ONLY"




2. How to Monitor Firmware update by USB Memory stick

[Description]


- Customer Condition

1. If the FW of USB must include a higher version more than version of current monitor, the updating feature will be normally worked.
2. User can choice whether or not performing a updating feature by the selection of user.



- Factory Condition (Aging On)

1. After searching the highest version in LG_MNT folder of USB, the updating feature has been automatically performed regardless of the current version of monitor.



[Common] When the updating feature is completed, monitor will turn off/on automatically after 3 seconds.

2. How to Monitor Firmware update by USB Memory stick

[Troubleshooting]

Q. The power turns off suddenly during updating FW in the status of progress bar, monitor does not work.

→ The power cable should be removed and reconnected after inserting an USB memory stick including FW in the determined folder.

[Note] After searching the highest version in LG_MNT folder of USB, the recovering feature has been automatically performed


→ After 40 seconds, the monitor will be normally recovered by itself.

Q. The power turns off suddenly during updating PD FW under 50% status of progress bar, how can I do?

→ After turning on the power, the updating will be progress again from the beginning percentage.

3. How to update DP Splitter Firmware

1. Connect the United jig to a PC and monitor





2. Enter the service menu using the following method

- 1) 'Setting > General > Information' and Operate the control key in the following order.
- 2) ▲ Press 2 times
- 3) ▼ Press 2 times
- 4) ► Press 5 times

3. Move to DP Split Update menu and turn OFF to ON.



4. Run the Gprobe program to update the firmware.



<fig-1>



3. How to update DP Splitter Firmware

Default settings-1





3. How to update DP Splitter Firmware

Default settings-2



3. How to update DP Splitter Firmware

Default settings-3




3. How to update DP Splitter Firmware

Setting commands in a text file

1. open Pegasus_ISP_UART.txt
2. Modify the red box in the figure to the right based on the firmware update file.
3. Save the file.



Pegasus_ISP_UART.txt



If the firmware is up to date,



Enter the filename here

3. How to update DP Splitter Firmware

If the connection is fine


1. enter 1 and press 'Execute' to see 0x01=0x80.



```
DP Splitter
0x01=0x80
Command Successful.
SetBuffer: Command Successful.
Reset: Command Successful.
Delay: Command Successful.
RAMWrite: Command Successful.
Run: Command Successful.
SetBuffer: Command Successful.
Delay: Command Successful.
Flash Manufacturer: 0xC2 ID: 0x13
SetDelay: Command Successful.
Erasing FLASH... Done.
```

To update the firmware,

1. navigate to the location of the text file and press 'Execute'.
2. The update will start, the screen will cycle off and on, and it will be completed when the following text appears.



```
batch D:\0000 ISP UART\Pegasus ISP UART.txt
PwrExtBlk-OFF
AppPwrOn
CHIP Power ON!!!
PwrExtBlk-ON
DP1.2 Pegasus Lib Version 0.25
19:42:12 Feb 25 2016
HDCP Repeater Lib Version 0.8
11:48:37 Sep 14 2015
HDMITX Lib Version 0.3
19:42:23 Feb 25 2016
MODEL: LG Medical Monitor
->PM:Pwr-On
```

4. How to update SDI Firmware

1. Please prepare an USB memory stick. (Except for the external HDD disk)

[Note] The file system must use the following format. (FAT32)


2. Please make a new folder like as LG_MNT of USB memory stick. (Ex. E:\₩LG_MNT)

3. Please copy SDI Software (2files :: FPGA_vx.xx_0xAAAA_yymmdd.bin , FW_vx.xx_0xAAAA_yymmdd.hex) to the folder of "LG_MNT".

[Note] The downloading feature does not work properly when using the download file from the company security.

4. Please insert the USB memory stick to the "SERVICE ONLY"

※ SDI firmware updates only support upgrading from a lower version to a higher version.




4. How to update SDI Firmware

5. Enter the Service mode


- 1) 'Setting > General > Information' and Operate the control key in the following order.
- 2) ▲ Press 2 times
- 3) ▼ Press 2 times
- 4) ► Press 5 times

6. In the USB Service Port menu, click the ► key



4. How to update SDI Firmware

7. SDI SW update pop-up message is shown



8. After showing FPGA Back up, FPGA, MCU Back up, MCU update, "Update process completed." message is shown about 3 sec . That means you are success to update

